

DR. SCOTT BULLOCH, DDS, MS
ORAL & MAXILLOFACIAL SURGEON

DR. BRANDON BULLOCH, DDS, MS
SURGICAL PROSTHODONTIST



PRACTICING SINCE 1994
PRECISION
DENTAL SPECIALTIES
ORAL SURGERY • PROSTHODONTICS

REFERRAL *Guide*

435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:
WWW.PRECISIONSTG.COM/REFERRINGDOCTORS



RESPECTED COLLEAGUES,

Over the years, we have had the privilege of working together on a multitude of cases, and each time, your referrals have reinforced the importance of a strong, united dental community. We firmly believe that a team approach is paramount in delivering comprehensive care that meets and exceeds our patients' expectations.

This informational book has been carefully designed to serve as a reference guide for everything you might need to know after referring a patient to us. It includes detailed information about our practice, the services we offer, and what your patients can expect during their treatment journey.

Our goal is to ensure that the transition from your care to ours is as seamless as possible. This guide provides insights into the procedures we perform, the technologies we employ, and the post-operative care guidelines that can assist you in coordinating with your patients. We believe that an informed and engaged patient is more likely to achieve the best possible outcomes.

In the coming weeks, we will be sending a copy of this guide to your office, and we encourage you to share it with your team and patients. Our hope is that this resource will empower you with the information you need and serve as a reference tool to address any questions or concerns that may arise when you refer patients to us.

Once again, we would like to extend our sincerest thanks for your ongoing support, trust, and collaboration. It is truly a pleasure to work alongside you, and we look forward to continuing to serve our community together.

If you have any questions or would like to discuss specific cases, please do not hesitate to reach out to us. We are always here to assist you in any way we can.

Thank you for being an integral part of our dental community, and for entrusting us with your patients' care.

Warm regards,

Dr. Scott Bulloch, DDS, MS
Oral & Maxillofacial Surgeon

Dr. Brandon Bulloch, DDS, MS
Surgical Prosthodontist

PRACTICING SINCE 1994
PRECISION
DENTAL SPECIALTIES
ORAL SURGERY • PROSTHODONTICS

TABLE OF CONTENTS

A LETTER FROM THE DOCTORS

PAGE 1

IMPLANT SYSTEMS

CURRENT SYSTEMS:

NOBEL PAGE 4

SOUTHERN PROMAX PAGE 10

STRAUMANN PAGE 13

PAST SYSTEMS:

SOUTHERN TRIMAX (RECENT PAST) PAGE 18

KEYSTONE TILOBEMAXX (RECENT PAST) PAGE 22

ASTRA PAGE 25

ZIMMER PAGE 29

FOR YOUR REFERENCE

DRIVERS FOR EACH SYSTEM PAGE 33

WARRANTY INFORMATION PAGE 34

LOCAL LABS REFERRING OFFICES USE PAGE 34

435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:

WWW.PRECISIONSTG.COM/REFERRINGDOCTORS

PRACTICING SINCE 1994

PRECISION

DENTAL SPECIALTIES

ORAL SURGERY • PROSTHODONTICS

CURRENT IMPLANT LINES



Nobel Active
Nobel Parallel
N1



ProMax



PAST IMPLANT LINES



Tilobemaxx



Tri-MAX



435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:

WWW.PRECISIONSTG.COM/REFERRINGDOCTORS



CURRENT IMPLANT *Lines*



NOBEL BIOCARE

THINGS TO KNOW:

The Nobel System is Color Coded:

Pink



Narrow Platform

Diameter: 3.75 mm

OR

N1 Diameter: 3.5

Yellow



Regular Platform

Diameter: 4.3 mm

OR

Diameter: 5.0 mm

OR

N1 Diameter: 4.0 mm

Blue



Wide Platform

Diameter: 5.5 mm

NON-NOBEL PERMANENT RESTORATIONS VOID PATIENT'S NOBEL REPLACEMENT WARRANTY.

435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:

WWW.PRECISIONSTG.COM/REFERRINGDOCTORS



PRACTICING SINCE 1994
PRECISION
 DENTAL SPECIALTIES
 ORAL SURGERY • PROSTHODONTICS

NOBEL BIOCARE IMPRESSION PARTS

Impression Copings & Analogs

Direct to Implant

Nobel Active/
Nobel Parallel

Analogs:

Internal conical connection
NobelActive, NobelParallel and
NobelReplace Conical Connection



NP	36697
RP	36698
WP	37879

Impression Copings:

Closed Tray	3.0	NP	RP	WP
Height	13.0 mm	13.0 mm	9.0 mm	13.0 mm
Ø 3.3 mm	36801			
Ø 3.6 mm		36538	36541	36540
Ø 3.8 mm	36803			
Ø 5.0 mm		36539	36543	36542 37851 37850
Ø 6.0 mm			36545	36544
Ø 6.5 mm				37853 37852

Nobel N1

Implant Replica Nobel Biocare N1™ TCC

NP	300962
RP	300963



Impression Coping Closed Tray Nobel Biocare N1™ TCC

	9 mm	13 mm
NP Ø 4.0	-	300957
RP Ø 4.5	300959	300958



NON-NOBEL PERMANENT RESTORATIONS VOID PATIENT'S NOBEL REPLACEMENT WARRANTY.

435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:

WWW.PRECISIONSTG.COM/REFERRINGDOCTORS



PRACTICING SINCE 1994

PRECISION

DENTAL SPECIALTIES

ORAL SURGERY • PROSTHODONTICS

NOBEL IMPRESSION PARTS

Impression Copings & Analogs

Direct to Base

Nobel N1



**IOS Implant Replica
Nobel Biocare N1™ TCC**

NP 300964

RP 300965




**Impression Coping Closed Tray
Nobel Biocare N1™ TCC**

	9 mm	13 mm
NP Ø 4.0	–	300957
RP Ø 4.5	300959	300958

Nobel On1

Analogs:



38731

38732

38733

Impression Copings:



38725

38726

38727

NON-NOBEL PERMANENT RESTORATIONS VOID PATIENT'S NOBEL REPLACEMENT WARRANTY.

435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:








WWW.PRECISIONSTG.COM/REFERRINGDOCTORS



PRACTICING SINCE 1994
PRECISION
 DENTAL SPECIALTIES
 ORAL SURGERY • PROSTHODONTICS



NOBEL DIGITAL PROCESS

Active & Parallel Scan Bodies

Implant connection	Platform	For single-unit implant restorations		For NobelProcera Implant Bridges		For models produced with IOS data**	
		Desktop scan body	IOS scan body	Desktop scan body for implant level	Desktop scan body for Multi-unit Abutment level	IOS model implant analog	IOS analog insertion tool
Internal conical connection NobelActive, NobelParallel and NobelReplace Conical Connection 	3.0	 2B-D	 IO 2B-D SA	 -	 -	 38188	 38200
	NP	2B-A	IO 2B-A SA	2B-A	2C-A	38189	38201
	RP	2B-B	IO 2B-B SA	2B-B	2C-A	38190	38203
	WP	2B-C	IO 2B-C SA	2B-C	2C-A	38191	38204

These engaging ELIOS Scan Bodies are What our Office Uses the Most!

N1 Scan
Bodies

Hand O-Mini		
IOS Healing Abutment Nobel Biocare N1™ Base		
NP	300991	
RP	300992	

On1 concept – componentry NP RP WP			
On1 IOS Healing Cap incl. handle			
4 mm			38739
4.5 mm	38735	38737	
5 mm			38738
6 mm	38734	38736	

On1 Scan
Bodies

If we know You Will be Scanning, We will Usually Place one of these at the time of Surgery.

435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:

WWW.PRECISIONSTG.COM/REFERRINGDOCTORS



PRACTICING SINCE 1994




PRECISION

DENTAL SPECIALTIES




ORAL SURGERY • PROSTHODONTICS

NOBEL ABUTMENTS

Active & Parallel: Single Abutments

			3.0			NP			RP			WP	
	Esthetic Abutment ¹	Height	1.5 mm	3.0 mm	4.5 mm	1.5 mm	3.0 mm	4.5 mm	1.5 mm	3.0 mm	4.5 mm	6.0×7.0 mm	7.0×8.0 mm
			36782	36783	36814	36665	36666	36249	36669	36671	36251	37827	37828
	Esthetic Abutment 15° ¹		36784	36785	36815	36667	36668	36250	36672	36673	36252		
	Universal Base ²	Height				1.5 mm	3.0 mm		1.5 mm	3.0 mm		1.5 mm	3.0 mm
						38213	38216		38214	38217		38215	38218

Multi-Unit Abutments:

		NP	RP	WP
	Multi-unit Abutment Plus ³			
	Height 1.5 mm	38878	38879	38880
	Height 2.5 mm	38881	38882	38883
	Height 3.5 mm	38884	38885	38886
	Height 4.5 mm		38887	
	17° Multi-unit Abutment Plus ³			
	Height 2.5 mm	38888	38889	38890
	Height 3.5 mm	38891	38892	38893
	30° Multi-unit Abutment Plus ³			
	Height 3.5 mm	38894	38895	
	Height 4.5 mm	38896	38897	

N1

Single Abutments:



**Universal Abutment
Nobel Biocare N1™ TCC**

	1.5 mm	3 mm
NP	300946	300947
RP	300948	300949



**Esthetic Abutment
Nobel Biocare N1™ TCC**

		0°	15°
NP	1.75 mm	300938	300939
	3 mm	300940	300941
RP	1.75 mm	300942	300943
	3 mm	300944	300945

N1

Multi-Unit Abutments:



**Multi-unit Abutment Xeal™
Nobel Biocare N1™ TCC**

	1.5 mm	2.5 mm	3.5 mm	4.5 mm
NP	301044	301045	301046	
RP	301047	301048	301049	301050



**17°/30° Multi-unit Abutment Xeal™
Nobel Biocare N1™ TCC**

	Angle	2.5 mm	3.5 mm	4.5 mm
NP	17°	301051	301052	
	17°	301053	301054	
RP	30°		301055	301056

NON-NOBEL PERMANENT RESTORATIONS VOID PATIENT'S NOBEL REPLACEMENT WARRANTY.



PRACTICING SINCE 1994

PRECISION

DENTAL SPECIALTIES

ORAL SURGERY • PROSTHODONTICS

NOBEL TOOLS

Required tools and torques

20 O-Mini 15 O-Mini

20/15 Ncm



Screwdriver Machine
OmniGrip™ Mini

20 mm	300852
28 mm	300853
36 mm	300854

Use this for N1 components linked to the base



Manual Torque Wrench Prosthetic	29165*
Manual Torque Wrench Adapter Prosthetic	29167

Our office's nickname for this driver is the "Star Wrench"

20 MUA

20 Ncm



Screwdriver Machine
Multi-unit

20 mm	29158
-------	-------

20 Uni 15 Uni

20/15 Ncm



Screwdriver Machine
UniGrip

20 mm	300852
28 mm	300853
36 mm	300854

Tighten to 20 Ncm

Screwdriver Nobel Biocare N1™

	28 mm	36 mm
Machine	300903	300904
Manual	301135	301136



On1 Base Screwdriver

Machine 21 mm	38752
---------------	-------

Tighten to 35 Ncm

NON-NOBEL PERMANENT RESTORATIONS VOID PATIENT'S NOBEL REPLACEMENT WARRANTY.

435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:

WWW.PRECISIONSTG.COM/REFERRINGDOCTORS



PRACTICING SINCE 1994

PRECISION

DENTAL SPECIALTIES

ORAL SURGERY • PROSTHODONTICS

SOUTHERN PROMAX



PROMAX6

Major Diameter:
6.0

Platform Diameter:
4.5



PROMAX7

Major Diameter:
7.0

Platform Diameter:
5.7



PROMAX8

Major Diameter:
8.0

Platform Diameter:
6.5



PROMAX9

Major Diameter:
9.0

Platform Diameter:
7.5

435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:

WWW.PRECISIONSTG.COM/REFERRINGDOCTORS






SOUTHERN PROMAX

Tools:

**THE PROMAX SYSTEM IS ACCESSED WITH A UNIVERSAL
 HEX DRIVER**



435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:


WWW.PRECISIONSTG.COM/REFERRINGDOCTORS



PRACTICING SINCE 1994
PRECISION
 DENTAL SPECIALTIES
 ORAL SURGERY • PROSTHODONTICS

SOUTHERN PROMAX

Impression Parts:



435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:

WWW.PRECISIONSTG.COM/REFERRINGDOCTORS






STRAUMANN

THINGS TO KNOW:

The Straumann System is Color Coded:



**Bone Level Tapered
SC**

Diameter:
2.9



**Bone Level Tapered
NC**

Diameter:
3.3

435-652-1445

office@precisionstg.com


REFER PATIENTS ONLINE:

WWW.PRECISIONSTG.COM/REFERRINGDOCTORS



PRACTICING SINCE 1994
PRECISION
DENTAL SPECIALTIES
ORAL SURGERY • PROSTHODONTICS

STRAUMANN



**Bone Level Tapered
SC**

Diameter:
2.9



**Bone Level Tapered
NC**

Diameter:
3.3

435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:
WWW.PRECISIONSTG.COM/REFERRINGDOCTORS






STRAUMANN

Tools:

THE PROMAX SYSTEM IS ACCESSED WITH A SCS SCREWDRIVER



435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:


WWW.PRECISIONSTG.COM/REFERRINGDOCTORS



PRACTICING SINCE 1994
PRECISION
 DENTAL SPECIALTIES
 ORAL SURGERY • PROSTHODONTICS

STRAUMANN

Impression Parts:



435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:

WWW.PRECISIONSTG.COM/REFERRINGDOCTORS





PAST IMPLANT *Lines*



SOUTHERN TRIMAX

THINGS TO KNOW:

The TRIMAX System is Color Coded:



TRIMAX 7

Major Diameter:
7.0
Platform Diameter:
5.7



TRIMAX 8 & 9

Major Diameter: 8.0	Major Diameter: 9.0
Platform Diameter: 6.5	Platform Diameter: 7.5

435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:

WWW.PRECISIONSTG.COM/REFERRINGDOCTORS



PRACTICING SINCE 1994

PRECISION

DENTAL SPECIALTIES

ORAL SURGERY • PROSTHODONTICS

SOUTHERN TRIMAX



TRI-MAX7-xx

Major Diameter:
7.0
Platform Diameter:
5.0



TRI-MAX8-xx

Major Diameter:
8.0
Platform Diameter:
6.0



TRI-MAX9-xx

Major Diameter:
9.0
Platform Diameter:
6.0

435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:
WWW.PRECISIONSTG.COM/REFERRINGDOCTORS



PRACTICING SINCE 1994
PRECISION
 DENTAL SPECIALTIES
 ORAL SURGERY • PROSTHODONTICS

SOUTHERN TRIMAX

Tools:

THE TRIMAX SYSTEM IS COMPATIBLE WITH THE NOBEL UNIGRIP

20/15 Ncm



**Screwdriver
Machine UniGrip**

20 mm	300852
28 mm	300853
36 mm	300854



Manual Torque Wrench Prosthetic	29165*
Manual Torque Wrench Adapter Prosthetic	29167

Our office's nickname for this driver
 is the "Star Wrench"

435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:


WWW.PRECISIONSTG.COM/REFERRINGDOCTORS



PRACTICING SINCE 1994
PRECISION
 DENTAL SPECIALTIES
 ORAL SURGERY • PROSTHODONTICS

SOUTHERN TRIMAX

Impression Parts:



435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:

WWW.PRECISIONSTG.COM/REFERRINGDOCTORS




PRACTICING SINCE 1994

PRECISION


DENTAL SPECIALTIES

ORAL SURGERY • PROSTHODONTICS


KEYSTONE TILOBEMAXX



Major Diameter:
7.0
Platform Diameter:
5.5



Major Diameter:
8.0
Platform Diameter:
6.5



Major Diameter:
9.0
Platform Diameter:
7.5

435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:
WWW.PRECISIONSTG.COM/REFERRINGDOCTORS






KEYSTONE TILOBEMAXX

Tools:

THE PROMAX SYSTEM IS ACCESSED WITH A QUAD DRIVER



435-652-1445

office@precisionstg.com


REFER PATIENTS ONLINE:

WWW.PRECISIONSTG.COM/REFERRINGDOCTORS



KEYSTONE TILOBEMAXX

Impression Parts:



435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:

WWW.PRECISIONSTG.COM/REFERRINGDOCTORS





ASTRA

THINGS TO KNOW:

The Astra Tech System is Color Coded:



EV 3.6 S

Diameter:
3.6



EV 4.2 S

Diameter:
4.2

435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:

WWW.PRECISIONSTG.COM/REFERRINGDOCTORS



PRACTICING SINCE 1994
PRECISION
DENTAL SPECIALTIES
ORAL SURGERY • PROSTHODONTICS

ASTRA



EV 3.6 S

Diameter:
3.6



EV 4.2 S

Diameter:
4.2

435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:
WWW.PRECISIONSTG.COM/REFERRINGDOCTORS






ASTRA

Tools:

**THE PROMAX SYSTEM IS ACCESSED WITH A UNIVERSAL
 HEX DRIVER**



435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:

WWW.PRECISIONSTG.COM/REFERRINGDOCTORS



PRACTICING SINCE 1994

PRECISION


DENTAL SPECIALTIES

ORAL SURGERY • PROSTHODONTICS


ASTRA

Impression Parts:

3.6




4.2



	3.6	4.2
Ø mm	3.6	4.2
Height mm	11	11
Order No.	25535	25537


	3.6	4.2
Height mm	16.5	16.5
Order No.	25544	25545

Truabutment Scan Body
REF#: AE36-SB



Before Scanning, Make
Sure that the Lab you
Use have Access to/Have
Downloaded the
Libraries to Support these
Scan Bodies OR Ask
What Scan Bodies they
Prefer you Use.

Truabutment Scan Body
REF#: AE42-SB



Traditional
Impression

OR

Scan

435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:

WWW.PRECISIONSTG.COM/REFERRINGDOCTORS








ZIMMER

THINGS TO KNOW:

The Zimmer Implant System is Color Coded:

Implant Diameter	Implant Platform	Color-coding
3.7mmD	3.5mmD	 Green
4.1mmD	3.5mmD	
4.7mmD	4.5mmD	 Purple
6.0mmD	5.7mmD	 Yellow

435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:

WWW.PRECISIONSTG.COM/REFERRINGDOCTORS








PRACTICING SINCE 1994

PRECISION

DENTAL SPECIALTIES

ORAL SURGERY • PROSTHODONTICS

ZIMMER

 Tapered Screw-Vent	Implant Diameter	Implant Platform
	3.7mmD	 3.5mmD
	4.1mmD	 3.5mmD
	4.7mmD	 4.5mmD
	6.0mmD	 5.7mmD

435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:
WWW.PRECISIONSTG.COM/REFERRINGDOCTORS






ZIMMER

Tools:

THE PROMAX SYSTEM IS ACCESSED WITH A UNIVERSAL HEX DRIVER



435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:

WWW.PRECISIONSTG.COM/REFERRINGDOCTORS



PRACTICING SINCE 1994
PRECISION
 DENTAL SPECIALTIES
 ORAL SURGERY • PROSTHODONTICS

ZIMMER

Impression Parts:


Traditional
Impression

OR

Scan



Implant Platform	Emergence Profile			
	3.5mmD	4.5mmD	5.5mmD	6.5mmD
3.5mmD	HLT3/3	HLT3/4	HLT3/5	•
4.5mmD	•	HLT4/4	HLT4/5	HLT4/6
5.7mmD	•	•	•	HLT5/6
3mm Extension Screw for Indirect Transfer	HLTE	HLTE	HLTE	HLTE
Replacement Retaining Screw	HLTS2	HLTS2	HLTS2	HLTS2
Implant Platform	Catalog No.			
	3.5mmD	4.5mmD	5.5mmD	6.5mmD
3.5mmD	IA3			
4.5mmD	IA4			
5.7mmD	IA5			



Truabutment
Scan Body
REF#: ZV3-SB,
ZV4-SB,
ZV5-SB

Before Scanning, Make Sure that the Lab you Use have Access to/Have Downloaded the Libraries to Support these Scan Bodies OR Ask What Scan Bodies they Prefer you Use.

435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:

WWW.PRECISIONSTG.COM/REFERRINGDOCTORS



Dr. Scott Bulloch, DDS, MS
Oral & Maxillofacial Surgeon



Dr. Brandon Bulloch, DDS, MS
Surgical Prosthodontist

DRIVER/ TORQUING FLOW CHART



N1

NP- 3.5 **PINK**
RP- 4.0 **YELLOW**

Implant Level OR Base

ONLY SCREW RETAINED TEMP CROWNS

Abutment Driver (Straight) 20 ncm



PROSTHETICS

Omnigrip M Driver
(Straight or ASC*) 20 ncm



Special Cases:

Multi-unit Prosthetic Screws- 15 ncm
Angled Multi-unit Abutments- 15 ncm
Straight Multi-Unit Abutments- 35 ncm



*Angled Screw Channel

Active/Parallel

NP- 3.5 **PINK**
RP- 4.5/5.0 **YELLOW**
WP- 5.5 **BLUE**

Implant Level Abutment

CEMENTED OR SCREW RETAINED
TEMP CROWNS

On1 Base Level

SCREW RETAINED TEMP CROWNS

Base Driver- 35 ncm



PROSTHETICS

Unigrip Driver (Straight) 35 ncm



Omnigrip Driver (ASC*) 35 ncm



FOR CEMENTED TEMPS:

Base Driver- Finger Tight



TRIMAX (Trilobe)

BLUE TRIMAX7- 5.7
GREEN TRIMAX8- 6.5
GREEN TRIMAX9- 7.5


Nobel Unigrip Driver
35 ncm



PROMAX (Hex)

PROMAX6- 4.5
PROMAX7- 5.7
PROMAX8- 6.5
PROMAX9-7.5

0.50 Hex Driver
35 ncm



TISSUE LEVEL

NNC- 3.3
RN- 3.3,4.1,4.5
WN- 4.8

Straumann Driver
35 ncm



BONE LEVEL

MOST COMMON

SC- 2.9
NC 3.3
RC 4.1,4.8


Straumann Driver
35 ncm



Phased Out Systems:


ASTRA

0.50 Hex
Driver



ZIMMER

0.50 Hex
Driver



KEYSTONE

TRIMAX7- 5.0
TRIMAX8/9- 6.0
Keystone Driver 35 ncm



Dr. Scott Bulloch, DDS, MS
Oral & Maxillofacial Surgeon



Dr. Brandon Bulloch, DDS, MS
Surgical Prosthodontist

WARRANTY

We have had great success with Nobel, Southern & Keystone's implant replacement warranties.

The few times we have had to replace implants through those companies, we have been successful at receiving replacements from companies.

The few times we haven't been successful receiving replacement implants have been with Nobel & Straumann when off-brand restoration brands have been used.

Our mutual patients pay a premium to use these premium implants. We recommend talking to the patients about their preference for branded restorations in order to maintain their implant warranty.

435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:

WWW.PRECISIONSTG.COM/REFERRINGDOCTORS

Dr. Scott Bulloch, DDS, MS
Oral & Maxillofacial Surgeon



Dr. Brandon Bulloch, DDS, MS
Surgical Prosthodontist

VETTED LOCAL LABS

***Beehive Dental Lab**

Joel- 435-632-0058

THIS IS OUR AO4 PARTNER LAB

***Color Country**

Devin- 435-628-1281

Set up to complete On1 & N1 Restorations

***Innovative Dental Studio**

Sean- 702-580-7900

Set up to complete On1 & N1 Restorations

***Digital Concepts**

(Previously Named Riverside Lab)

435-628-1444

Set up to complete On1 & N1 Restorations

Epic Dental Lab

Scott- 801-756-1117

Set up to complete On1 & N1 Restorations

***Local to St. George Area**

435-652-1445

office@precisionstg.com

REFER PATIENTS ONLINE:

WWW.PRECISIONSTG.COM/REFERRINGDOCTORS